

ITEM QTY	DESCRIPTION	REF
1	3 GLOVE ISOLATOR SHELL	ST.STL
1	FILTER SUPPORT FRAME	ST.STL
1	FAN CASSETTE CHAMBER	ST.STL
1	FAN HOUSING	ST.STL
2	EBM DOUBLE INLET FAN	FAN-422
1	180 x 600 VHP INLET CHAMBER	ST.STL
1	180 x 600 VHP FUMIGATION PANEL	ST.STL
1	1365 x 533 x 66 INLET HEPA FILTER	H14
2	INTERNAL LIGHT LENS	SILICONE
2	LIGHT LENS SEAL/LOCKING STRIP	SILICONE
4	ADJUSTABLE FEET	ST.STL
1	RRI MAGNETIC HATCH	ST.STL
1	CLARUS PORT	ST.STL
1	PASS BOX VHP FILTER CLAMP FRAME	ST.STL
1	15SER HATCH VHP FUMIGATION PANEL	ST.STL
1	1460 x 320 x 124 PB EXT FILTER	H14
1	1460 x 320 x 71 PB INT FILTER	H14
18	SET CORZAN VHP PIPEWORK	PBW-012/08
1	RRI POLYPROPYLENE SLIDING TRAY	ST.STL
21	CONTROL PANEL	ST.STL
8	VISOR CLAMPS	ST.STEEL
1	VISOR FRAME	COATED ALUM
1	INTERLOCK CONTROL FOOTSWITCH	ODD-028
24	INTERLOCK CONTROL FOOTSWITCH	ODD-028
1	WASTE PORT LID	ST.STL
1	WASTE PORT CHUTE	ST.STL
1	VISOR SCREEN	ST.STL
1	LOWER FRAME	CLEAR ACRYLIC
1	MEASURE FRAME	ST.STL
1	UPPER FRAME	ST.STL
2	BLEED FILTER	H14
2	BLEED FILTER CLAMP	ST.STL
2	DOP/PARTICLE PORTS	ST.STL
1	HANGING RAIL C/W 12 HOOPS	FTC-923
36	SET LAMINAR FLOW SCREENS	CLEAR ACRYLIC
1	PRE-FILTER	ST.STL

REV	DATE	ISSUED FOR INFORMATION ONLY	DESCRIPTION
REV 1	01/12/09	ISSUED FOR INFORMATION ONLY	

ENVAIR Ltd
York Avenue, Huddersfield,
England, BS4 4HX
Tel: +44(0)1708 208416
Fax: +44(0)1708 831957
Email: Technical@envair.co.uk
Internet: www.envair.co.uk

Client: ENVAIR GmbH

Contract No. 28714 Works Order

Drawn: S.THEOBOLD Date: 30/11/09

Scale: 1:10 Check:

Drawing No. Issue

ISO1-230 1 OF 1

